

---

---

**NEISSERIA GONORRHEA BY PCR**

For the detection of *Neisseria gonorrhoea* (GC) by Roche Amplicor PCR nucleic acid amplification method in endocervical, male urine, and male urethral specimens. For other sites refer to Culture, GC.

**Specimen Requirements:****Endocervix and male urethra:**

Complete instructions, collection supplies, and required transporters included in each kit available from St. Marys Laboratory Services.

**Male urine:**

Patient must not have urinated during the previous 2 hours. Collect **first void urine**, a minimum of 10 mL up to 50 mL in a sterile, leak proof container. Refrigerate until pickup.

**Reference Range:**

Not applicable

**Report Time:**

24 - 72 hours

**Days(s) Set Up:**

Monday through Friday

**CPT Codes(s):**

GC by PCR 87591